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Indications

Instillamed® is a medical device for the temporary replacement of the GAG layer on the bladder epithelium in:

- all forms of abacterial chronic and chronic recurrent cystitis
- chronic recurrent cystitis
- interstitial cystitis
- radiogenic cystitis

Description

The inner wall of the bladder is lined with an impermeable protective layer, the so-called GAG (glycosaminoglycan) layer. This GAG layer is made up primarily of chondroitin sulphate and hyaluronic acid. Damage to the GAG layer compromises its protective function and its ability to shield off irritating substances such as potassium ions, microcrystals, bacteria and proteins, which have a harmful effect on the bladder epithelium.

Instillamed® has been specially developed for the temporary replacement of the protective GAG layer on the bladder epithelium and contains highly purified chondroitin sulphate, manufactured using a special process, and hyaluronic acid.

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Instillamed® Instructions for use. Please read carefully.

Composition

50 ml Instillamed® contains 800 mg sodium hyaluronate, 1000 mg chondroitin sulphate.

Posology and method of administration

The bladder should first be emptied completely. The entire contents of the Instillamed® pre-filled syringe is then instilled into the bladder via the catheter. For this purpose, the syringe can be attached to an appropriate catheter using the Luer lock connector. If a catheter without a Luer lock connector is being used, the syringe can be attached to the catheter using the adapter supplied. The solution should remain in the bladder for as long as possible, however for at least 30 minutes. At the start of treatment, depending on the course of symptoms, one instillation per week for 4 weeks of treatment is recommended. The treatment should then be given for up to 6 months, at monthly intervals.

Contraindications

Instillamed® should not be used in children under the age of 12 years of age or during pregnancy, as there is no experience in the use of Instillamed® in these patient groups. Do not use Instillamed® if you have macrohaematuria or if you have suffered an injury to the urethra and/or the bladder.

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Juli 2017

Notes

Do not use Instillamed® if you have a known allergy to one of its constituents or if you are allergic to fish protein. Do not use if the pre-filled syringe is damaged. Intended for single use only.

As in any urological use, a urinary tract infection may occur due to the insertion of the catheter.

This preparation must be used only by a medical practioner who has been specially trained in the treatment.

Storage

Do not store above 25 °C. Do not freeze. Keep Instillamed® out of the reach and sight of children.

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FARCO-PHARMA GmbH, Gereonsmühlengasse 1-11, 50670 Köln, Germany, Tel: +49(0)221/594061
Melisana AG, Grüngasse 19, CH-8004 Zürich, Tel.:+41 44 247 72 00
ELVIM Ltd., Kurzemes prosp. 3, Office 513, Riga, LV-1067, Latvia, Tel.: +371/6139744
Medeco B.V., Alexander Flemingstraat 2, 3261 MA Oud Beijerland, Nederlande, Tel.: +31/186634400