

I. Video

II. FAQ


Purenum GmbH | Fahrenheitstraße 1 | D-28359 Bremen
Tel. +49 421 2208 170 | www.purenum.com | vigilanz@purenum.com

STERILE

5°C / 41°F

30°C / 86°F

MD

en **mediNiK® - Basic** I. Training video II. FAQ

Two-component hydrogel to remove kidney stone fragments

The following equipment is needed to use mediNiK®: a sterile syringe containing 5 mL of the „K1“ blue component (basic substance) and a sterile syringe containing 11 mL of the „K2“ yellow component (activator). Each sales unit contains five syringes each of K1 and K2.

Intended Use: mediNiK® is used to enclose kidney stones and concretions in a hydrogel for removal during ureterorenoscopy (URS). The use of mediNiK® allows fragments to be removed from the kidney that otherwise could not be removed with a grasping instrument due to their small size.

Users: For use by healthcare professionals only.

Contraindications: Hypersensitivity or allergic reaction to any of the ingredients. The contraindications and restrictions of ureterorenoscopy apply (coagulation disorders, anticoagulation therapy, sepsis, untreated urinary tract infections, pregnancy, kidney disease, intolerance to endoscopic procedures, marked hyperplasia of the prostate, ureteral stenosis).

Pregnancy and lactation, children and adolescents: There are no data available regarding the use of mediNiK® in pregnant or breast-feeding women or in children and adolescents.

Combination with accessories: There are no known restrictions on the choice of endoscopes, grasping instruments, and ureteral access sheaths (UAS).

Before use of mediNiK®:

A. Make sure that there are no contrast agents, blood, or other substances in the kidney that could affect visualisation and/or hydrogel formation.

B. Before removing the syringe cap, release the pressure point by pressing the plunger down. This ensures that the syringe can be easily and evenly emptied.

C. For best results, the kidney stones should be gathered at the deepest point of a renal calyx.

During use of mediNiK®:

1. **Steer** the endoscope until the kidney stone fragments are visualised. Flush with sterile saline solution (0.9% NaCl) to ensure good visualisation.

2. Stop flushing. Administer the K1 blue basic substance slowly by pulsing it through the working channel of the endoscope until the fragments to be removed are completely enclosed in the K1 blue substance. If more than one renal calyx is affected, repeat this step for each calyx.

Tipp: Only use as much of the K1 blue basic substance as you need. This will be transformed into the hydrogel following the administration of the second component, the K2 yellow activator.

3. Flush the working channel of the endoscope carefully with 0.9% NaCl.

Caution! Avoid turbulence/dilution of the K1 blue substance by flushing the working channel too vigorously. Slightly withdraw the endoscope.

4. Administer the K2 yellow substance carefully and slowly through the now-flushed working channel. Move the endoscope slowly towards the affected renal calyx/calyses.

Caution! Avoid turbulence the K1 blue substance by too rapid administration. Regardless of how much K1 blue basic substance was used, the entire contents of the K2 yellow activator syringe (11 mL) must be administered.

5. Wait for approx. 3 minutes until the hydrogel has sufficiently set. Waiting an inadequate time, or diluting the K2 yellow activator through flushing, will result in a softer hydrogel.

6. Grasp and remove the hydrogel clot carefully using the grasping instrument. The area can now be flushed again, this time more gently.

Caution! The grasping instrument can break the hydrogel if it is held too tightly.

After use: Make sure that all hydrogel clots were removed. If any hydrogel clots do remain, these are – usually – dissolved as part of diuresis. The hydrogel can also be actively dissolved via administration of a chelating agent (e.g. EDTA).

Repeating the procedure: If the procedure is repeated, the entire quantity of K2 yellow activator used previously must be flushed from the kidney.

Single use: mediNiK® is intended for single use and must not be resterilised.

Disposal: Syringes containing the K1 blue and K2 yellow substances can be disposed of in accordance with the lowest level of disposal standards at the hospital. The hydrogel must be disposed of separately as it has come into contact with mucus and/or urine.

Note: All serious events occurring in relation to the product must be reported to the manufacturer and the competent authorities of the Member State in which the user and/or patient is located (vigilanz@purenum.com).

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