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Instructions for use. Please read carefully before use.

Endosgel® Chlorhexidine Digluconate 0.05%

Sterile lubricant with physiological conductivity and ancillary disinfectant properties.

Endosgel® is used in adults, adolescents (12 - 18 years) and children (2 - 11 years).

Contains: Chlorhexidine digluconate 0.05%, methyl-4-hydroxybenzoate, propyl-4-hydroxybenzoate, sodium-S-lactate, hydroxyethyl cellulose, purified water.

Intended purpose/indication for use: **Endosgel®** is a sterile lubricant with ancillary disinfectant properties to facilitate the insertion of catheters and instruments into the urethra and bladder, particularly during endourological procedures, electro-coagulation and resection and for intermittent (self-) catheterisation. **Endosgel®** is used in diagnostic investigations of disorders of the urethra and bladder, for instance cystoscopy, in electrocauterisation of bleedings in the urethra or bladder, and to facilitate catheterisation in bladder voiding disorders.

Properties: **Endosgel®** leaves instrument lenses with a clear view and has physiological conductivity (≈ 7 to 10 mS/cm). The positively charged chlorhexidine digluconate has an ancillary disinfectant action by penetrating the cell wall of microorganisms. It thereby alters their surface structure, which leads to reorientation of the cell wall molecules and to a disturbance of osmotic equilibrium. This ultimately results in destruction of the microorganism's cell wall. **Endosgel®** therefore possesses ancillary disinfectant properties with regard to urothelial pathogens, whereby urinary tract infections can be avoided.

Interactions: The action of chlorhexidine can be impaired by protein, blood, soap and other anionic substance. Strong bases and strong acids interact with chlorhexidine; precipitates occur at a pH > 8 whereas hydrolysis takes place at a pH < 3. Chlorhexidine interacts with strong oxidising agents, so do not use **Endosgel®** at the same time with strong oxidising agents on the same areas.

Side effects: Skin manifestations can occur in the event of hypersensitivity; very rarely (<1/10 000 persons), contact dermatitis, bladder irritation or allergic reactions up to anaphylactic shock to chlorhexidine digluconate can occur.

The physician will decide on the required treatment measures.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (contact details of competent authorities via link: <https://cesportal.hma.eu/>).

Contraindications: **Endosgel®** must not be used in patients with hypersensitivity/allergy to alkyl-4-hydroxybenzoates (parabens), chlorhexidine digluconate or another ingredient.

Warnings and precautions:

Chlorhexidine has a negative influence on the microflora.

When **Endosgel®** is used, the benefit-risk relationship must be considered by a physician, particularly

- in patients with a weakened immune system, as secondary infections can occur.
- when injury or changes of the mucosa are suspected.
- when local infection of the body orifice is suspected.
- when **Endosgel®** is used in children, adolescents and elderly patients.

Endosgel® should be used by medical professionals except in the case of intermittent (self-) catheterisation. In this case, the user must be instructed by the physician in handling and dosage. In the case of uncertainty, the treating physician should be contacted directly.

Protect the product from overheating during use.


There is no adequate experience or studies of the safety of use of **Endosgel®** in pregnancy and lactation period.

Use with other products: There is no evidence of undesirable interactions or impaired function between **Endosgel®** and instruments and medical devices commonly used in the context to the intended purpose of **Endosgel®**.

The choice of a suitable catheter should be made by a physician. If a lubricant is required also with coated catheters, you should discuss the use of **Endosgel®** with your doctor.

It is recommended to clarify the use of **Endosgel®** together with active medical devices with the respective manufacturer.

Dosage instructions:

Preparation: Remove the syringe from the blister. Before removing the sealing cap, a possible pressure point should be loosened by lightly pressing the plunger (press the sealing cap against firm resistance). This enables easy and uniform emptying of the syringe. The sealing cap is then removed, as far as possible without touching the syringe tip.

Please note that there is an air bubble in the syringe. This does not influence Endosgel functioning.

Note: For the gel to reach its full effectiveness, there should be a waiting period of 5-10 minutes after instillation.

Use: Unless prescribed/specified otherwise:

a) Transurethral insertion of catheters and instruments: Instill the contents of the 6ml or 11ml syringe. After the usual cleaning and disinfection of the surrounding skin and/or mucosal areas, e.g., the external genitalia and external urethral orifice, **Endosgel®** is instilled slowly into the urethra. The entire urethra including the external sphincter must be coated with a lubricant film for gentle introduction of catheters and instruments. The volume depends on the anatomy and/or the type of procedure. The single-use syringe containing 6ml **Endosgel®** is suitable for use in women and

the 11ml one in men. When an instrument is introduced for the first time, use of another syringe of Endosgel can be indicated clinically, e.g., because of the anatomy (in male patients, for example, when there is the condition after known urethral strictures, especially in case of enlarged middle prostate lobe or when there is the condition after Vïa Falsa; in female patients, for example, suffering from urethra-hypospadias or adipose patients, in which cases the urethra is difficult to identify) or because of the type of procedure (e.g., electroresection and coagulation). If an instrument or catheter is changed during a procedure, it is advisable based on clinical experience to instil a further syringe of **Endosgel®**. The physician will decide on the need for further instillation during a procedure. For all uses, not more than 7 syringes with 11ml units or 10 syringes with 6 ml units should be used per day.

b) 2-4-ml from the 6 ml syringe are recommended in children (2 – 11 years). The volume depends on the age and/ or the type of procedure and on the anatomy. There are no special recommendations on dosage, and the dosage is determined by the physician. For dosing of **Endosgel®**, please remove the air bubble prior to instillation by holding the opened syringe upwards vertically and, as soon as the air bubble is at the upper end of the syringe, expelling it by lightly pressing the plunger. A small amount of gel is also expelled. There are markings on the syringe plunger. These act as an approximate 1ml guide during instillation of smaller volume units. The graduation does not have any measurement function.

Endosgel® must not be used in children under 2 years.

c) Intermittent (self-) catheterisation:

Endosgel® should be used by medical professionals except in users who have to carry out intermittent (self-) catheterisation. In this case, the user must be instructed by the physician in handling and dosage. In the case of questions or uncertainty, the patient should consult a physician. Make hygienic preparations (wash and disinfect your hands, ideally wear gloves, prepare an area where germs are reduced as much as possible) and place all necessary aids ready to hand.

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Instructions on (self-) catheterisation can be sent on request. Additional information on (self-) catheterisation can also be found on homepage of FARCO-PHARMA GmbH (www.farco.de) in the section for patients under Catheterisation in German and English.

The choice of a suitable catheter should be made by a physician.

During use you should not eat, drink or smoke. **Endosgel®** must not be introduced into other body orifices than stated in the instructions for use (e.g., mouth, eye, ears or nose). Otherwise, wash out immediately with clear water.

Additional information:

Notes: **Endosgel®** is intended for single use and must not be resterilised. Discard leftover gel and do not use the empty syringe for other purposes as germ transmission can occur.

Storage: The syringe should be stored in the original packaging (blister in carton) at room temperature, dry and protected from light (do not exceed a storage temperature of 25°C, do not freeze). The product is sterile as long as the blister is sealed and undamaged. Do not use the syringe if the blister is opened or damaged.

Shelf life: **Endosgel®** must not be used after the stated expiry date. The expiry date refers to the last day of the stated month.

Disposal: **Endosgel®** can be disposed of with household waste according to national guidelines, but should preferably be handed into a hazardous waste collection point. Do not dispose of expired or unused product in the toilet or wash basin. Products must not get into waste water or sewerage. If disposed of with household waste, ensure that there can be no improper access to this waste. If you are unsure, ask your physician or pharmacist how to dispose of the product.

Contents: Packs of 10 and 100 single-use syringes containing 6/11 ml.

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Medeco B.V.
Melisiana AG
FARCO-PHARMA GmbH, Gereonsmühlengasse 1-11, 50670 Köln; Telefon: +49(0)221/594061;