
1.Language

mediNiK®

I. Video
I. Training video
II. FAQ
II. FAQ

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STERILE

B

Two-component hydrogel for the removal of kidney stone fragments

The following components are required for one application of mediNiK®: One sterile syringe containing 5 mL of the blue component K1 (base substance) and one sterile syringe containing 11 mL of the yellow component K2 (activator). One sales unit contains five K1 syringes and five K2 syringes.

Intended purpose: mediNiK® is used to enclose kidney stones and concretions in a hydrogel for removal during ureteronoscopy (URS). The use of mediNiK® allows fragments to be removed from the kidney that otherwise could not be removed with a grasping instrument due to their small size.

CE region: To expand the therapeutic benefit, component K1 and/or K2 can be mixed with water-soluble substances/active ingredients prior to application.

Intended users: For professional use only.

Indication: Urolithiasis

Contraindications: Hypersensitivity or allergic reaction to any of the ingredients.

Composition: Blue Dextran, calcium chloride, disodium hydrogen phosphate, glycerol, potassium dihydrogen phosphate, riboflavin, sodium alginate, sodium chloride, water.

Pregnancy, breastfeeding, children and adolescents: There is no experience with the use of mediNiK® in pregnant or breastfeeding women, children or adolescents.

Combination with accessories: There are no known restrictions regarding the selection of endoscopes, stone retrieval baskets, flexible suction sheaths (FANS) and ureteral access sheaths (UAS).

Before using mediNiK®:

A. Ensure that there is no contrast medium, blood or other substances in the kidney that could affect visibility and/or hydrogel formation.

B. Before removing the syringe cap, depress the plunger rod slightly to overcome the initial plunger resistance. This ensures smooth and even emptying of the syringes.

C. For an optimal result, the kidney stones should be collected at the lowest point of a renal calyx.

Use of mediNiK® with a stone retrieval basket:

1. Navigate the flexible or semirigid ureteroscope under direct vision to the kidney stone fragments. Irrigate with sterile physiological saline (0.9% NaCl) to ensure good visibility.

2. Stop irrigation! Slowly apply the blue base substance K1 in pulses through the working channel of the endoscope, until the fragments to be removed are completely surrounded by blue K1. If more than one renal calyx is affected, repeat this step separately for each calyx.

Tip: Use only as much of the blue base substance K1 as necessary. It subsequently forms the hydrogel upon addition of the second component, the yellow activator K2.

3. After application of K1, a separating layer of NaCl must be created between the components to prevent gel formation in the working channel. To do so, withdraw the endoscope to the entrance of the calyceal system and slowly introduce 1-2 mL of NaCl into the working channel.

Caution! Avoid dispersing/diluting the blue K1 by irrigating the working channel too vigorously. Withdraw the endoscope slightly.

4. Carefully and slowly apply the yellow K2 through the working channel. Slowly retract the endoscope towards the affected renal calyx/calyces.

Caution! Avoid dispersing the blue K1 by adding K2 too quickly. Regardless of the amount of blue base substance K1 used, the entire contents of the yellow activator K2 syringe (11 mL) should be applied.

5. Insert the stone retrieval basket (e.g. a Dormia basket or a tipless basket) through the working channel of the endoscope.

6. Carefully grasp and remove the hydrogel clot with the stone retrieval basket. Irrigation may now be resumed at low flow.

Caution! The basket may cut through the hydrogel if it is closed too tightly.

Use of mediNiK® with a suction sheath: See FAQ QR link and/or direct QR link.

After use: Ensure that all hydrogel clots have been removed. If any hydrogel clots remain, they will generally dissolve as a result of diuresis. The hydrogel can also be actively dissolved by adding a chelating agent (e.g. EDTA).

Repeating the procedure: If the procedure is repeated, all previously used yellow K2 must be completely flushed out of the kidney.

Single use: mediNiK® is intended for single use only and must not be resterilized.

Transport: 5°C – 40°C

Disposal: The syringes containing blue K1 and yellow K2 may be disposed of in accordance with the lowest applicable hospital-internal disposal standard. The hydrogel must be disposed of separately because it has come into contact with mucous membranes and/or urine.

CE region: If water-soluble substances/active ingredients have been added, disposal must be carried out in accordance with the requirements applicable to those substances/active ingredients.

Note: Any serious incident occurring in relation to the device must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established (vigilanz@purenum.com)

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